AMILON

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AUTOMATIC STARCHMETER FOR CONTROLLING THE RIPENESS OF APPLES AND PEARS

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EN The **AMILON** is an automatic starchmeter that photographs an equatorial section of fruits that have been previously exposed to an iodine-potassium iodide solution. Innovative software processes the image and expresses a precise measurement of the starch stage. This instrument is particularly suitable for all cases in which it is necessary to determine the index of ripeness in apples and pears with a quick measurement that is accurate and immune from erroneous interpretation.

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ADVANTAGES

Objective and automated measurement of starch hydrolysis (starch stage), a universally used assay as an index of ripeness in apples and pears; reproducible instrumental measurement immune from erroneous interpretation of the visual reading; quick measurement on single piece and/or sample of max 12 pieces; measurement feasible for all variants without having to refer to specific figures.

OPERATION PRINCIPLE

The starch stage is highlighted through reaction with an iodine-potassium iodide solution. The equatorial section of the fruit thus treated is photographed, then the image is processed and the starch value is computed and displayed, for each single piece and as an average over the sample (max 12 pieces). The starch value can be displayed in accordance with commonly used reference scales of 1÷10 (internationally recommended), 1÷5, 1÷8 or 1÷1000.

**SPECIFICATION**

Size	400x450x600 mm
Weight	6 Kg
Minimal system requirements for operating software	1.8 GHz processor or higher MS Windows 95/98/2000/Me/XP or NT 4.0 256 MB RAM memory or higher VGA screen card with 8 MB memory or higher 150 MB of available disk space USB 2.0 port or higher, preferably USB 2.0